Willson Gelatin Sponge
Sterile and haemostatic gelatin sponge
Instructions for use
Manufacturer:

WILLPHARMA
- SINCE 1924 -

Beechavenue, 6
1119 PT Schiphol-Rijk
The Netherlands

Phone: +31 (0)20 497 65 51
E-mail: info@willpharma.com
www.wilsson.com

Instructions for Use – Medical device

What is Wilsson?

Wilsson Gelatin Sponge is a sterile gelatin sponge from porcine origin, insoluble in water, easy to handle, re-absorbable and which can be cut to the desired size without crumbling. Its purpose is to facilitate the stopping of bleeding of venous or arterial capillaries, and in cases where other conventional haemostatic procedures (stitching, compression, binding, etc.) cannot be carried out. The sponge consists of 99.3 % porcine gelatin and 0.7% of sodium lauryl sulphate and is whitish with uniform porosity. Wilsson Gelatin Sponge is sterilized by gamma radiation.

How does Wilsson work?

Once saturated with blood, the product promotes coagulation, on average within 3 minutes, by serving as a coadjuvant for controlling local haemorrhage. If properly utilized and in suitable amounts, Wilsson Gelatin Sponge is reabsorbed within 3-4 weeks. The amount of product to be used and its reabsorption depends on the extension of the bleeding surface and the type of tissue on which it is applied.

Presentation

Wilsson exists in 4 sizes:

- Wilsson Gelatin Sponge Standard (ref. 20925NL-20): individually packaged in a double package consisting of an inner and an outer pouch. Dimensions*: 80 x 50 x 10 mm – Box of 20 units
- Wilsson Gelatin Sponge Spécial (ref. 20945NL-20): individually packaged in a double package consisting of an inner and an outer pouch. Dimensions*: 80 x 50 x 1 mm – Box of 20 units
- Wilsson Gelatin Sponge Dental (ref. 20966NL-50): individually packaged in a blister pack. Dimensions*: 10 x 10 x 10 mm – Box of 50 units
- Wilsson Gelatin Sponge Tampon (ref. 21005NL-16): individually packaged in a blister pack. Dimensions*: 80 x 30 0 mm – Box of 16 units

*Sizes may vary with +/- 2 mm.

Indications

Wilsson Gelatin Sponge can be used dry or soaked in sterile saline solution, for surgical interventions, except for ophthalmic ones, through simple application or by exerting slight pressure. Wilsson Gelatin Sponge can be cut to the desired size without crumbling, using sterile scissors, and it has a low adhesiveness with respect to forceps and surgery instruments.

Wilsson Gelatin Sponge Standard can be used in parenchymal, thoracic, ENT (ear, nose throat) surgery and neurosurgery ^(1,3).

Wilsson Gelatin Sponge Spécial is a thinner type of sponge for use in vascular surgery to reduce bleeding when reconstructing blood vessels ⁽²⁾. Wilsson Tampon is indicated for use in anal and rectal surgery ⁽²⁾. Wilsson Gelatin Sponge Dental can be used in dentistry following tooth extraction and prior to filling cavities resulting from the removal of mandibular or maxillary cysts ⁽⁴⁾.

Contraindications

- Do not use Wilsson Gelatin Sponge for closing skin cuts as it could interfere with healing.
- Do not use Wilsson Gelatin Sponge in intravascular cavities due to the risk of embolization.
- Do not use Wilsson Gelatin Sponge for embolization procedures.
- Do not use Wilsson Gelatin Sponge in patients with associated allergy towards the product.
- Do not use Wilsson Gelatin Sponge on infected wounds.
- Do not use Wilsson Gelatin Sponge in cases of strong arterial haemorrhage, in cases of blood stagnation or stagnation of other organic liquids where the source of haemorrhage is entirely occluded.
- Do not use Wilsson Gelatin Sponge in post-partum interventions and for controlling such haemorrhage.

Beware

- Wilsson Gelatin Sponge is a medical device retaining parts of porcine animal tissue. The physician should inform the patient of this.
- Wilsson Gelatin Sponge is supplied sterile and disposable. Unused parts of the product should be discarded. Any opened and unutilized packaging should be discarded in biological waste.
- The use of Wilsson Gelatin Sponge does not authorize abandonment of normal haemostatic intervention.
- The safety and efficiency of Wilsson Gelatin Sponge in ophthalmic interventions is not verified.
- The safety and efficiency of Wilsson Gelatin Sponge on children and expectant women is not verified.

How to use Wilsson?

Carefully read the present instructions for use before using the product. Ensure that the packaging is whole before using. Sterility is guaranteed only if the packaging is not damaged or wet.

Use the product by the expiry date indicated on the packaging. Always be cautious when it comes to aseptis whenever handling the product. Always use the minimal amount of product to obtain haemostasis.

Wilsson Gelatin Sponge can be used dry as is or soaked with saline solution. Apply the product on the area subject to haemorrhage by applying slight pressure until haemostasis is obtained. When haemostasis is not reached in 3 minutes due to a larger bleeding, you can apply a second Wilsson Gelatin Sponge.

After obtaining haemostasis, if possible, remove the product from the site of application to avoid the risk of displacing the product and possible compression of the adjacent structures. If the product cannot be removed upon obtaining haemostasis, carefully eliminate the surplus used product and place it in biological waste.

Remove Wilsson Gelatin Sponge from the site of application when used in bone holes, on the spinal cord and/or on the nerve or optical chiasm or in the area surrounding or adjacent thereto.

Remove Wilsson Gelatin Sponge from the site of application if used in neurosurgery or in the areas surrounding or adjacent to the nerves given that by increasing the volume the product may cause compression leading to damage of nerve tissue.

Precaution measures - Warnings

- When in use in closed or narrow tissue spaces preliminary compression of the product is advisable to prevent its increase volume-wise, due to contact with organic fluids, causing damage.
- When indicated to fill a cavity to obtain haemostasis, Wilsson Gelatin Sponge should not be used this way unless all surplus product is removed after achieving haemostasis
- Use Wilsson Gelatin Sponge in minimal amounts to reach haemostasis and always remove any surplus product from the site of application. When surplus material is not in contact with blood, reabsorption can take longer, which might require a reintervention.
- Do not use Wilsson Gelatin Sponge together with autologous blood recovery systems, given that it has been proven that particles of haemostatic agents could pass through the transfusion filters of the purification systems.
- Do not use Wilsson Gelatin Sponge together with methyl methacrylate cement. It has been proven that the porous structure weakens the cement used to fix the prosthesis to the bone surface.
- There have been cases of incomplete reabsorption and loss of hearing due to the use of haemostatic sponges in ear-drum-plastic operations.
- Do not use Wilsson Gelatin Sponge in cases of coagulation disorders.
- The association of Wilsson Gelatin Sponge with pharmacologically active substances or other substances shall be evaluated by the surgeon who takes responsibility thereof.
- The safety and efficiency of Wilsson Gelatin Sponge in urology surgical procedures are not verified. Wilsson Gelatin Sponge should not be left in the renal pelvis, in renal calyces, in the urethra or in the ureter to avoid the possible formation of calculi.
- Wilsson Gelatin Sponge is not recommended for treating cases of epistaxis given that the product could be inhaled thus causing respiratory problems.
- Wilsson Gelatin Sponge is for single use. Do not re-use material due to loss of sterility.
- Wilsson Gelatin Sponge Sterile Gelatin Sponge should not be re-sterilized by heat, because heating may change absorption time. Ethylene oxide is not recommended for re-sterilization because it may be trapped in the interstices of the foam. Although, not reported for Wilsson Gelatin Sponge, the gas is toxic to tissue, and in trace amounts may cause burns or irritation.

Storage

Store Wilsson Gelatin Sponge in a dry place between 15°C and 25°C. Keep Wilsson Gelatin Sponge away from direct sunlight.

Symbols used on the product label

Consult instructions For use	Single use	Do not re-sterilize
Caution	Do not use if package is damaged	Sterilized using irradiation
Use by date of Expiry Date (yy-mm-dd)	Batch code or lot number	Date of manufacture (yy-mm-dd)
Catalog number	Magnetic resonance safe	

References

- Cho KS et al., 2015 Comparative analysis of Cutanplast and Spongostan nasal packing after endoscopic sinus surgery: a prospective, randomized, multicenter study. Eur Arch Otorhinolaryngol. 2015;272:1699-705.
- Doria C, et al. Phase 3, randomized, double-blind study of plasma-derived human thrombin versus bovine thrombin in achieving haemostasis in patients undergoing surgery. Curr Med Res and Opin. 2008;24:785-94.
- Liang J, et al. Using tranexamic acid soaked absorbable gelatin sponge following complex posterior lumbar spine surgery: A randomized control trial. Clin Neurol Neurosurg. 2016;147:110-114.
- Napenas JJ, et al. Review of postoperative bleeding risk in dental patients on antiplatelet therapy. Oral Surg Oral Med Oral Pathol Oral Radiol. 2013;115:491-9.
- Kobatake K, et al. Effect on haemostasis of an absorbable haemostatic gelatin sponge after transrectal prostate needle biopsy. Int Braz J Urol: official journal of the Brazilian Society of Urology. 2015;41:337-43.

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Références

- Cho KS et al., 2015 Comparative analysis of Cutanplast and Spongostan nasal packing after endoscopic sinus surgery: a pro-spective, randomized, multicenter study. Eur Arch Otorhinolaryngol. 2015;272:1699-705.
- Doria C, et al. Phase 3, randomized, double-blind study of plasma-derived human thrombin versus bovine thrombin in achieving haemostasis in patients undergoing surgery. Curr Med Res and Opin. 2008;24:785-94.
- Liang J, et al. Using tranexamic acid soaked absorbable gelatin sponge following complex posterior lumbar spine surgery: A randomized control trial. Clin Neurol Neurosurg. 2016;147:110-114.
- Napenas JJ, et al. Review of postoperative bleeding risk in dental patients on antiplatelet therapy. Oral Surg Oral Med Oral Pathol Oral Radiol. 2013;115:491-9.
- Kobatake K, et al. Effect on haemostasis of an absorbable haemostatic gelatin sponge after transrectal prostate needle biopsy. Int Braz J Urol: official journal of the Brazilian Society of Urology. 2015;41:337-43.

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