Hersteller

lege artis Pharma GmbH + Co. KG
Postfach 60, D-72132 Dettenhausen
Breitwasenring 1, D-72135 Dettenhausen
Telefon +49 (0) 71 57 / 56 45 - 0
Fax +49 (0) 71 57 / 56 45 50
E-Mail: info@legeartis.de
www.legeartis.de

Bezeichnung

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(GB) Expert information
CALCINASE EDTA-Solution

CE 0124

Manufacturer

lege artis Pharma GmbH + Co. KG
P. O. Box 60, D-72132 Dettenhausen
Breitwasenring 1, D-72135 Dettenhausen
Telephone +49 (0) 71 57 / 56 45 - 0
Fax +49 (0) 71 57 / 56 45 50
E-Mail: info@legeartis.de
www.legeartis.de

Name

CALCINASE EDTA-Solution
For rinsing of root canals
For dental use only

Composition

Purified water, disodium edetat and sodium hydroxide

Intended purpose

CALCINASE EDTA-Solution works by chemically complexing calcium ions and is used for rinsing and widening the root canal and removing smear layer during preparation.

Indication

For root canal treatment - for rinsing and widening the root canal and for removing smear layer during the preparation.

Other applications

In addition, it can be used for finding root canals and for removal of harmful calcium deposits on sensitive instruments.

Contra-indications

Allergy to disodium edetate
Wide open apical foramen

Side-effects

Side-effects are not known for intended use. Very long contact times or rinsing with large amounts of EDTA-Solution may result in softening of the root dentin and increased dentin permeability. Note: please also see

Mode of application

If the substance passes through the apex, irritations of the periapical tissue are possible.

In the event of serious incidents in connection with the product, please contact the manufacturer and the responsible authority immediately.

Notices and precautions

Avoid swallowing. After swallowing, rinse mouth.

May cause skin and eye irritation. In case of skin contact, wash with plenty of water and soap. If eye contact occurs, rinse with water for a few minutes, remove any existing contact lenses if possible. In case of persistent irritation, seek medical advice/medical attention.

Independent of the mode of application in the root canal, suitable instruments and a qualified method (slow application without pressure, removing the rinsing solution by suction, protection of gingiva and oral mucosa by use of rubber dam) have to be used during the whole application. The resulting smear layer (dissolved dental substance, residual pulp and microorganisms) is rinsed out of the root canal after each change of instrument size. For this, alternating rinsing with sodium hypochlorite solution is recommended. The rinsing intensity has to be adapted in accordance with the amount of used CALCINASE EDTA-Solution. Finally, the root canal will be rinsed with an inert solution (e.g. physiological saline solution).

The safety data sheet can be downloaded from www.legeartis.de or requested by E-Mail from sicherheitsdatenblaetter@legeartis.de.

Mode of application

Removal of smear layer

To remove the smear layer developed during preparation, a large-volume rinsing of at least 1 minute in the prepared root canal is necessary.

An alternating rinsing with sodium hypochlorite solution is to be recommended. Note: Rinse protocols in German are available, ask info@legeartis.de

Widening of root canals

If the root canals are too narrow, the root canal must be rinsed with CALCINASE EDTA-Solution by means of e.g. a pipette or a syringe or a cotton wool thread soaked with CALCINASE EDTA-Solution must be introduced and then the root canal must be prepared with the one size larger root-canal file. Rinsing with CALCINASE EDTA-Solution and the further preparation must be repeated until the root canal is sufficiently widened.

Note: Rinse protocols in German are available, ask info@legeartis.de

Finding of root canals

In case of difficult-to-find root canal entrances CALCINASE EDTA-Solution is to be instilled by means of e.g. a pipette or a syringe into the cavities of the lower teeth and to be introduced by means of a cotton wool pellet into the upper teeth. The canal entrances become clearly visible within a few minutes.

Removal of harmful calcium deposits on sensitive instruments

In case of a calcified Airotor water nozzle CALCINASE EDTA-Solution is filled by means of e.g. a pipette or a syringe directly into the emptied water nozzle until the liquid passes out of the nozzle of the angular head. After 2 -3 minutes rinse with water.

In obstinate cases, it is also possible to help with the thin wire supplied with each turbine from the side of the nozzle orifice.

Sieves and sensitive instruments, decalcify best by brushing with or short soaking in CALCINASE EDTA-Solution.

Note: Please follow the care instructions of the instrument manufacturer!

Further notice

CALCINASE EDTA-Solution contains 20% of disodium edetate. Disodium edetate is used for the acid-free dissolution of dentine for the preparation of the root canal and serves the removal of the smear layer, especially in combination with sodium hypochlorite solution.

One molecule of disodium edetate forms a slightly water-soluble, relatively stable chelate complex with a calcium ion. As a result, apatite is dissolved, hard tooth substance is demineralised, and the existing smear layer is removed. The demineralising effect is time-dependent, in wide (large-volume or widely prepared) root canals known, in narrow root canals and in the apical third part of the root it is expected to be significantly less effective. Disodium edetate increases the dentine permeability.

The effect of disodium edetate is self-limiting, i.e. the demineralisation only lasts until all molecules of disodium edetate have bound to their calcium ion. The chelate complex can migrate into the dentine and via the root canal into the periapical tissue, traces remain in the dentine. Consequently, care should be taken to ensure sufficient rinsing during and at the end of the preparation to prevent too much dissolution of apatite, demineralization and softening of tooth substance.

Rinsing with chlorhexidine solutions immediately after using CALCINASE EDTA-Solution produces a white precipitate. Therefore, it is necessary to either rinse a little longer or rinse with saline solution between.

Do not pour back residues of the solution into the bottle.

Keep out of reach of children.

Shelf life

CALCINASE EDTA-Solution shall not be used after the expiry date.

Administrative form and package sizes

50 ml Solution Item number 0032301

200 ml Solution Item number 0032332

500 ml Solution Item number 0032331

Date of revision
2024-08

Symbol of "Expiry date"

Symbol of "Batch number"

Symbol of "Pay attention to the directions for use"

Symbol of "Manufacturer"

Symbol of "Catalogue number"
Shows Item number of the manufacturer


Symbol of "Medical Device"

Symbol of "Unique device identification"

Symbol of "upper temperature limit"


1. Handling of the ESD-syringe filling system with Luer or Luer Lock syringe

Remove the cap

Connect the syringe

Withdraw the desired volume

Remove the syringe

Close the cap

2. Pouring out the solution (without syringe) is also possible.

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